# Euroline

## Eurokat

### Eurokat – A Valuable Alternative for Food Analysis

Eurokat high performance polymer phases were especially developed for the separation of organic acids, carbohydrates, alcohols and even complex mixtures of these compounds. Eurokat is a sulfonated cross-linked styrene-divinylbenzene copolymer available in several ionic forms (H, Ca, Pb). This particular cation exchanger is characterized by an outstanding density of functional groups, making it the ideal choice for your ion exclusion and ligand exchange chromatography.

Typical particle size:	10 µm
Particle distribution:	± 1 µm
Cross-linkage:	about 8%
Max. pressure:	up to 300 bar

#### Stability

Eurokat polymer columns are extremely stable over the whole pH range. This is one striking advantage compared with silica-based phases where lifetime, especially in the higher pH range, is limited. Most importantly, the Eurokat phases show extraordinary mechanical stability and are not affected by aqueous solvents.

#### Application Areas

Three different ionic species are available. Eurokat H is especially designed for determinations of organic acids and even complex mixtures of acids, sugars and alcohols as well as sugar alcohols. Eurokat Ca and Eurokat Pb are the phases of choice for determinations of smaller carbohydrates (DP < 4). The application area of Eurokat also includes the determination of biochemical species such as sugars, artificial sweeteners or acidic metabolic and degradation products. These compounds often play a key role in the analysis of food and beverage products and in clinical applications including:

- Carbohydrates, acids in fruit and soft drinks
- Sugar substitutes in diet beverages
- Detection of illegal additives in wine
- Sweets and candies
- Dairy products (aging, detection of butyric acid)
- Food preservatives
- Krebs cycle acid intermediates
- Urine (uric acid, hippuric acid)
- Monitoring of fermentation processes

#### Choice of Eluents

Outstanding separation results are achieved isocratically with fresh double-distilled water for Pb or Ca type Eurokat columns, while Eurokat H columns perform best when using inorganic acids. Such aqueous eluents produce no harmful waste as with organic solvents such as acetonitrile which necessitate expensive disposal.

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## Choice of Chromatographic Conditions


The Eurokat stationary phase is designed to be applicable to a wide range of diverse chromatographic conditions. All Eurokat columns can be used at temperatures up to 90°C under diverse pH conditions (e.g. using 0.001 – 0.1 N H<sub>2</sub>SO<sub>4</sub>). For example, the best separation of sugars can be achieved using an acidic aqueous mobile phase at a temperature between 60 to 90 degrees Celsius.

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
### Separation of Monosaccharides

1. saccharose
2. glucose
3. fructose
4. mannose
5. xylose
6. sorbitol

Column: Eurokat Pb, 10 µm, 300 x 8 mm ID  
 Eluent: H<sub>2</sub>O  
 Flow rate: 0.6 ml/min  
 Temperature: 80°C  
 Detection: RI

## KNAUER Column Hardware

We design and manufacture HPLC column hardware ranging from 2 mm ID to 62 mm ID under strict quality control. A specially treated inner surface ensures consistent column packing and high column stability. A wide range of column lengths from 5 mm up to 300 mm are available. An easily exchangeable integrated precolumn is available upon request.

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Dr. Ing. Herbert Knauer GmbH